

Combination cold and compression therapy to treat and prevent injury.

put a wrap on pain



VITALWEAR

RECOVERYWRAP®  
EQUINE SYSTEM

## A breakthrough in injury recovery and prevention. An effective, affordable and veterinarian-approved treatment that's easy enough for everyday use.

Cumbersome ice packs and hosing were once the only option to fight pain following injury or surgery. The RecoveryWrap Equine System combines the benefits of ice with constant compression therapy, delivered through specially designed wraps that fit easily onto a horse's body. Used in 30-minute cycles at regular intervals, RecoveryWrap Equine eases suffering as it accelerates the healing of damaged tissue. Applied after a rigorous workout, the RecoveryWrap Equine System can actually prevent strains from becoming full-blown injuries.

### **SIMPLE**

Designed for easy, everyday use

### **PORTABLE**

Rechargeable battery pack enables therapy when and where horses and trainers need it

### **FLEXIBLE**

Cold and compression therapies

### **IDEAL FOR**

Post-workout cool down; reducing inflammation & pain; recovery from tendon or ligament injury; post-surgical recovery; treating lameness; controlling muscle damage



### **Wide-ranging applications — One versatile solution**

Bowed tendon • Chronic lameness • Fetlock lameness • Preventative maintenance • Hock injury  
Post-op recovery • Suspensory demitis • Tendon sheath effusion • Many other conditions



VITALWEAR



put a wrap on pain

### RecoveryWrap Equine Bundles

01140 RecoveryWrap Equine System for Localized Leg Treatment (includes: RecoveryWrap Equine System, right & left front leg wraps, right & left rear leg wraps, right & left hock wraps, carrying bag, single & Y-splitter tubing set and other accessories)

01150 RecoveryWrap Equine System for Full Leg Area Treatment (includes: RecoveryWrap Equine System, right & left full front leg wraps, right & left full rear leg wraps, carrying bag, single & Y-splitter tubing set and other accessories)

01160 RecoveryWrap Equine System for Spine and Upper Leg Treatment (includes: RecoveryWrap Equine System, poll wrap, back wrap, sacroiliac wrap, stifle wrap, carrying bag, single & Y-splitter tubing set and other accessories)

### Systems

order no.

00400 RecoveryWrap System 15Vdc

### Wraps

00801 Right Front Leg

00802 Left Front Leg

00803 Right Rear Leg

00804 Left Rear Leg

00805 Right Hock

00806 Left Hock

01032 Right Full Front Leg

01033 Left Full Front Leg

01024 Right Full Rear Leg

01025 Left Full Rear Leg

01034 Poll

01035 Stifle

01036 Back

01037 Sacroiliac

00001 6" Universal Wrap

00118 8" Universal Wrap

### Accessories

00857 Equine Carrying Bag

00948 Tote Bag

00949 Long Elastic X 1

00950 Short Elastic X 2

00994 Y-Splitter tube set

00034 Single tubing set

00262 Power cord

00496 RW Wall Adapter Power Supply

00491 RW Battery Pack

### SYSTEM SPECIFICATIONS:

<b>PHYSICAL</b>	
size (approximate)	8" diameter x 11.8" height (20.3 cm x 30.0 cm)
weight (dry)	5 lbs (2.27 kg)
fluid tether	8' (2.44 m)
<b>CONTROL SYSTEM</b>	
type	manual user control
<b>THERMAL SYSTEM</b>	
range	40° - 70° F (4.4° - 21.1° C)
contrast therapy transition time	1 min typical from extremes in temperature range
operating period	2-6 hours typical continuous use
<b>CIRCULATING SYSTEM</b>	
reservoir capacity	0.8 gal (2.9 L)
reservoir fluid	ice water
flow rate thru VW	5.5 gph (21 Lph)
maximum pressure	18 psi
<b>ELECTRICAL SYSTEM:</b>	
voltage	15 VDC
power	15 W
leakage current	under 300 µA
<b>OPERATING ENVIRONMENT:</b>	
atmospheric pressure	525 - 795 mmHg
humidity	30% - 70% relative non-condensing
temperature	10° C - 40° C
<b>BATTERY PERFORMANCE:</b>	
80 minute charge	2 hours of use (approx.)
overnight charge	2 1/2 hours of use (approx.)
<b>TRANSPORT &amp; STORAGE ENVIRONMENT:</b>	
atmospheric pressure	179 TO 760 mmHg
humidity	10% - 95% relative non-condensing
temperature	0° C - 50° C

VitalWear, Inc.  
 384 Oyster Point Blvd., Suite 16  
 South San Francisco, CA 94080-1968  
 Tel: 1-800-553-4081  
 www.vitalwear-equine.com

For more information or to order VitalWear Equine Systems, please call 1-800-553-4081 or visit [www.vitalwear-equine.com](http://www.vitalwear-equine.com).

00835 Rev. A